PCT


WORLD INTELLECTUAL PROPERTY ORGANIZATION International Bureau



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification 6:	4.2	(11) International Publication Number: WO 00/02045
G01N 33/50	A2	(43) International Publication Date: 13 January 2000 (13.01.00)
(21) International Application Number: PCT/BES (22) International Filing Date: 6 July 1999 (CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL,
(30) Priority Data: 98870149.6 6 July 1998 (06.07.98) (71) Applicant (for all designated States except US): SCREEN S.A. [BE/BE]; Avenue des Becassines 7, Brussels (BE).	I EUR	Published Without international search report and to be republished upon receipt of that report.
 (72) Inventors; and (75) Inventors/Applicants (for US only): DUPRIEZ, [BE/BE]; Avenue de Roodebeek 229, B-1030 (BE). DETHEUX, Michel [BE/BE]; Rue Ch Loupoigne 13/102, B-1348 Louvain-La-Neuv PARMENTIER, Marc [BE/BE]; Chaussée d'Uc B-1630 Brussels (BE). (74) Agents: VAN MALDEREN, Eric et al.; Office Van M Place Reine Fabiola 6/1, B-1083 Brussels (BE). 	Brussonarles we (Bl cole 30	

(54) Title: HIGH-THROUGHPUT SCREENING DIAGNOSTIC AND/OR DOSAGE METHOD OF AN AGONIST AND/OR AN ANTAGONIST FOR A CALCIUM-COUPLED RECEPTOR



(57) Abstract

The present invention is related to a diagnostic and/or dosage method of an agonist and/or an antagonist and/or a modulator for a calcium—coupled receptor and/or channel and/or any other calcium—coupled protein, comprising the following successive steps: disposing the agonist and/or the antagonist and/or the modulator upon a solid support; incubating one or more cell(s) expressing apoaequorin and said calcium—coupled receptor with coelenterazine in order to reconstitute an active aequorin in said cell(s); adding to said solid support one or more of said cells; and obtaining the measurement of an emitted light by said cell(s). The present invention is also related to the diagnostic and/or dosage device intended for the method according to the invention.

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
ΑT	Austria	FR	France	LU	Luxembourg	SN	Senegal
ΑU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav	TM	Turkmenistan
BF	Burkina Faso	GR	Greece		Republic of Macedonia	TR	Turkey
BG	Bulgaria	HU	Hungary	ML	Mali	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MN	Mongolia	UA	Ukraine
BR	Brazil	IL	Israel	MR	Mauritania	UG	Uganda
BY	Belarus	IS	Iceland	MW	Malawi	US	United States of America
CA	Canada	IT	Italy	MX	Mexico	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NE	Niger	VN	Viet Nam
CG	Congo	KE	Kenya	NL	Netherlands	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NO	Norway	zw	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's	NZ	New Zealand		
CM	Cameroon		Republic of Korea	PL	Poland		
CN	China	KR	Republic of Korea	PT	Portugal		
CU	Cuba	KZ	Kazakstan	RO	Romania		
CZ	Czech Republic	LC	Saint Lucia	RU	Russian Federation		
DE	Germany	LI	Liechtenstein	SD	Sudan		
DK	Denmark	LK	Sri Lanka	SE	Sweden		
EE	Estonia	LR	Liberia	SG	Singapore		

5

10

HIGH-THROUGHPUT SCREENING DIAGNOSTIC AND/OR DOSAGE METHOD OF AN AGONIST AND/OR AN ANTAGONIST FOR A CALCIUM-COUPLED RECEPTOR

Field of the invention

The present invention is related to a high15 throughput screening diagnostic and/or dosage method and
device of an agonist and/or an antagonist for a calciumcoupled receptor and the agonist and/or antagonist of said
calcium-coupled receptor identified by said method and
device.

20

Background of the invention and state of the art

A lot of G-protein-coupled receptors (GPCR) trigger, upon binding of an agonist, a transient increase in intracellular calcium concentration. This variation acts as an internal secondary messenger and is an important modulator of many physiological mechanisms (reviewed by Rink (1990), Tsunoda (1993) and by Santella & Carafoli (1997)). Measurement of intracellular calcium concentration in cells expressing a GPCR can thus be used to monitor the efficacy of activation of a GPCR by various compounds known - or suspected - to be a ligand for this GPCR.

The activation of other receptors such as ions channels may also induce an intracellular calcium concentration.

2

Changes in calcium concentration can be detected by several means and methods, like the use of fluorescent dyes (for example: fura-2, fluo-3, fluo-4 and indo-1).

However, Ca++ sensitive dyes have limitations. Activation of the dyes with an excitation beam requires complicated and expensive instruments and limits the use of the plastic labware such as microtiter plates.

Another method for intracellular calcium 10 concentration measurement is the use of cell overexpressing a GPCR and apoaequorin, such as described by Sheu et al. (1993). In this system, cells expressing apoaequorin are incubated with coelenterazine, which is the co-factor aequorin. During this incubation, of enters the cell and conjugates 15 coelenterazine apoaequorin to form aequorin, which is the active form of the enzyme. Upon incubation of the cells with an agonist of the GPCR, intracellular calcium concentration increases. This increase leads to the activation of the catalytic 20 activity of aequorin, which oxidises coelenterazine and yields apoaequorin, coelenteramide, CO2 and light. Once the photon has been emitted, the complex must dissociate and apoaequorin must recombine with a new coelenterazine molecule to be able to emit light again. Thus, in this system, measurement of light emission following agonist 25 addition reflects its ability to activate the GPCR and thus to increase intracellular calcium concentration. Because light is emitted only during 20 to 30 seconds after activation of the GPCR, recording of the emitted light must 30 be performed during the few seconds following agonist addition to the cells. This flash-type signal is due to the fact that (1) intracellular calcium increase triggered by GPCR is only transient and (2) as mentioned earlier, after

3

oxidation of coelenterazine, apoaequorin must recombine with coelenterazine to be able to emit light again.

The Patent Application EP-0341477 teaches the jellyfish photoprotein aequorin expression of 5 mammalian cell system by cloning gene pAQ440 specifying the biosynthesis of the aequorin into an expression vector a mammalian cell system, subjecting the plasmid of resulting plasmid to transfection and producing the photoprotein aequorin in the mammalian cell.

Patent US-5,422,266 describes a encoding apoaequorin protein included in a vector capable of expressing the apoaequorin in E. coli.

10

20

The Patent US-5,714,666 describes mammalian cell lines or transgenic animals expressing apoaequorin and 15 a receptor involved in the modulation of intracellular calcium. This document also describes a method of measuring intracellular calcium comprising adding coelenterazine cofactors to said mammalian cells expressing apoaequorin and measuring photoemission where emission of photons is indicative of intracellular calcium concentration.

However, the methods of the state of the art require firstly the spreading of cells from a mammalian cell line expressing apoaequorin on a solid support (for example a 96-well plate), secondly the addition of the 25 coelenterazine cofactor upon the cells and incubation to reconstitute a functional aequorin, thirdly the preparation agent affecting a receptor involved in the modulation of intracellular calcium concentration, and its addition to the prepared cells, and finally the measurement 30 of the photoemission.

Furthermore, as mentioned above, emitted only during 20 to 30 seconds after activation of the GPCR. Therefore, the recording of the emitted light

4

must be performed during the few seconds following agonist addition to the cells.

Therefore, the methods used in the state of the art are not adequate for a detection based upon high-throughput screening level, which usually need luminometer(s) and require the use of microtiter plates for the testing of thousands of compounds.

Aims of the invention

The present invention aims to provide a method and means, which do not present the drawbacks of the state of the art, for detecting biologically active substances, especially agonists and/or antagonists for calcium-coupled receptors.

A main aim of the present invention is to provide such method and means which allow the detection of biological active substances at a high-throughput scale, which could be adapted to specific recipients such as microtiter plates without requiring the modification of the high-throughput screening device.

Another aim of the present invention is to provide an easy and non-expensive method that could be easily automated.

25 Summary of the invention

The present invention is related to a high-throughput screening diagnostic and/or dosage method of an agonist and/or an antagonist for a "calcium-coupled" receptor, comprising the following successive steps:

- 30 disposing an agonist and/or an antagonist (preferably of a molecule) of said receptor upon a solid support,
 - incubating one or more cell(s) expressing apoaequorin and said "calcium-coupled" receptor with coelenterazine

WO 00/02045

5

PCT/BE99/00087

in order to reconstitute an active aequorin by said cell(s),

- adding to said solid support one or more of said cells,
 and
- 5 obtaining the measurement of an emitted light by said cell(s).

The term "calcium-coupled" receptor means any receptor (such as a G-coupled receptor or an ion channel) whose activation (by an ion, a known or unknown agonist or antagonist molecule) may increase an intracellular calcium concentration in the cell comprising said receptor, preferably in its cytoplasmic membrane.

The terms "disposing... upon a solid support" means the step of putting said compound into a recipient such as a microtiter plate or any other solid support without requiring any (covalent or other) fixing of said compound to said solid support.

Advantageously, the solid support is a microtiter plate, preferably a 96-well microtiter plate.

- Advantageously, the cell expressing apoaequorin and the calcium-coupled receptor is a cell expressing a G-coupled receptor and possibly one or more protein(s) intended to ensure a coupling of the receptor to the calcium pathway.
- Preferably, said protein is selected from the group consisting of a natural $G\alpha 16$ or $G\alpha 15$ protein, a chimeric G-protein resulting from a fusion between two different G-proteins or a phospholipase $C\beta 2$ protein.

The measurement of the emitted light is advantageously obtained with one or several luminometer(s), possibly equipped with several dispensers and measurement heads.

6

The present invention is also related to a high-throughput screening diagnostic and/or dosage device intended for the diagnostic and/or dosage method according to the invention, said device comprising the following 5 elements:

- a recipient, preferably a microtiter plate, more preferably a 96-well microtiter plate,
- a medium containing cell(s) expressing apoaequorin and a calcium-coupled receptor,
- 10 a medium containing coelenterazine, and
 - means (such as one or more luminometer(s) equipped with one or several dispensers and measurement heads) for detecting and possibly quantifying an emitted light by said cell(s).
- Advantageously, the device according to the invention comprises means for the automatic performance of the successive steps of the diagnostic and/or dosage method according to the invention.

A last aspect of the invention is related to

the agonist and/or the antagonist of a calcium-coupled
receptor identified by the method or the device according
to the invention.

The present invention will be described in details in the following non-limiting examples, in reference to the enclosed figures.

Brief description of the drawings

Figure 1 shows a series curves representing the intensity of the emitted light by cells as a function of time for each well of a 96-well plate injected with cells expressing the CCR5 receptor, apoaequorin and Gal6. The scaling is the same for all the graphs. Recording of

WO 00/02045

5

7

the signal was performed for 30 seconds. Ligand concentrations are increasing from column 1 towards column 12. All measurements were performed in duplicate: lines A and B: ligand is RANTES; lines C and D: ligand is MIP-1 α ; lines E and F: ligand is MIP-1 β ; lines G and H: ligand is derivative A of RANTES.

PCT/BE99/00087

- Figure 2 represents the dose-response curve for different agonists of the CCR5 receptor which represent RLU (integration of emitted light on 30 seconds) according to the logarithm of the final concentration of the ligand.
- Figure 3 represents the dose-response curve for different agonists of the 5HT-2B receptor.
 - Figure 4 represents the dose-response curve of the light emitted from K562 cells expressing CCR3 and aequorin as a response to the activation of the receptor by eotaxin.
- 20 Figure 5 represents the dose-response curves for cells expressing aequorin and $G\alpha 16$ and (panel A) the orexin 1 receptor or (panel B) the orexin 2 receptor.
- Figure 6 represents the dose-response curve for different antagonists of the 5HT-2B receptor.

Description of a preferred embodiment of the present invention

Detection of agonistic activities by means of 30 mammalian cell lines expressing apoaequorin and a GPCR requires the measurement of the emitted light to be performed just after placing the cells in contact with the potential agonist. This can easily be measured at low

throughput using a single-tube luminometer. However, up to now, this biological system could not be used at a high-throughput scale. Indeed:

8

- the necessity to measure light just after placing the cells in contact with the agonist to be tested compels 5 use a luminometer equipped with a build-in dispenser. For example, due to the short duration of light emission, it is impossible to inject the drugs to be tested on the cells placed in the 96 wells while the luminometer 10 plate is outside the subsequently record emitted light with the plate into the luminometer. Even if the plate could rapidly (i.e. than 15 seconds) be placed into less luminometer after injection of the drugs to be tested, 15 current apparatus do not allow the measurement of light from the 96 wells before the extinction of the flash signal of aequorin, as these luminometers are not equipped with 96 detectors.
- 20 luminometers equipped with a build-in dispenser only
 20 allow to inject a single solution into the 96 wells,
 making it impossible to inject a different drug in
 each well. Moreover, the washing of the dispenser
 before each measurement, for the injection of another
 drug in the next well, is time-consuming and thus is
 25 not suitable for the high-throughput scale. The same
 problem occurs with devices equipped with 6 dispensers
 (e.g. the "Microbeta Jet" from EG&G Wallac) as the
 6 dispensers only deliver a single solution.

The present invention provides a method for performing high-throughput screening of drugs binding to GPCR by the use of mammalian cell lines expressing apoaequorin and a GPCR and by the use of a conventional luminometer. Following this method, the solutions to be tested for (ant)agonistic activities are placed in the

9

wells of a 96-well plate. Cells expressing apoaequorin and a GPCR are detached from the culture plate (or collected suspension cultures) and are incubated coelenterazine to reconstitute active aequorin. These are then maintained in suspension with a magnetic stirrer and the cell suspension is injected, well by well, on the solutions of potential agonist to be tested. Light emission is then recorded for 1 (alternatively up to 30 or more) seconds. This method, by injecting the same cell suspension in each of the 96 wells, avoids the need of washing the 10 dispenser between each measurement and allows to perform 96 measurements of agonist-induced aequorin light emission in 15 minutes or less with a single dispenser luminometer. Alternatively, it allows to perform 96 measurements of agonist-induced aequorin light emission in 2 minutes or 15 less with a luminometer equipped with 6 dispensers and measurement heads (e.g. with the "Microbeta Jet" from EG&G Wallac).

This method thus allows to perform highthroughput screening (10 000 samples/day) with mammalian cell lines expressing apoaequorin and a GPCR and by the use of a conventional luminometer. This reduces the screening time and the amount of drugs needed for each measurement.

This system also allows to perform a functional screening with very few (down to 5000 or less) cells per measurement.

The injection of the cells into the wells containing the agonists did not increase the background of the measurement (that could for example have originated from cell breakage, releasing aequorin molecules from the cells into the culture medium, where the calcium concentration would have triggered the emission of light from aequorin). A signal-to-noise ratio above 50 was commonly obtained with this system of cell injection.

10

The method according to the invention is suitable for performing high-throughput analysis of GPCR or other calcium-coupled-receptor stimulation by known or potential agonists by means of cells expressing 5 receptor and apoaequorin. These cells may express apoaequorin in the cytoplasm, as described by Sheu et al. (1993) or Button and Brownstein (1993) or may express apoaequorin in the mitochondria, by means of the addition of a mitochondrial targeting sequence to the aequorin, as 10 used by Stables et al. (1997) or in any other part of the cell. These cells may also express proteins intended to ensure coupling of the over-expressed receptor to the calcium pathway. These may be the natural $G\alpha 16$ or $G\alpha 15$ proteins (Milligan et al., 1996), chimeric G proteins resulting from a fusion between two different G proteins (Komatsuzaki et al., 1997), phospholipase $C-\beta 2$ (Park et al., 1992), or any other "universal coupling" protein. Once the cells have been prepared and loaded with coelenterazine, they can be used for several hours (at least 9 hours). The load in coelenterazine and the 20 intensity of the light emitted by the cells upon agonist stimulation lasts is stable for this period of time.

Examples

25 Example 1

A CHO cell line expressing the chemokine CCR-5 receptor, the Gα16 coupling protein and apoaequorin was established. Cells were cultivated as a monolayer in HAM'sF12 medium containing 10% Foetal bovine serum (FBS).

30 On the day of the experiment, the culture medium was removed and cells were incubated for 5 min at room temperature in PBS-EDTA (phosphate buffered saline solution without calcium, supplemented with 5 mM EDTA). Cells were

11

detached from the culture vessel by shaking the culture plate by hand and by pipetting up and down. Cells were centrifuged and the supernatant was removed to eliminate the EDTA; the pellet was resuspended in HAM's F12 culture medium without FBS and with 0.1% Bovine Serum Albumin. Cells were counted by means of a Thomas cell, were centrifuged again and were resuspended in HAM's F12 culture medium without FBS and with 0.1% Bovine Serum Albumin at a concentration of 5.10^6 cells/ml. Coelenterazine (or a derivative of it , e.g. Coelenterazines f, h, n, cp or hcp, from Molecular Probes Inc.) at 500 μ M in methanol was added to the cell suspension at a final concentration of 5 μ M. The cell suspension was then stored in the dark at room temperature for 3 to 5 h, with shaking every 15 to 30 min to maintain the cells in suspension.

Series of dilutions of known ligands were prepared in HAM's F12 culture medium without FBS and with 0.1% Bovine Serum Albumin and 50 μ l of each of these solutions were placed in the wells of a 96-well plate. The cell suspension was diluted 5 times with medium HAM'sF12 20 without FBS and with 0.1% Bovine Serum Albumin and was placed in a glass or plastic container protected from light by wrapping it with aluminium foil. A magnetic stirring bar was added to the suspension and a magnetic stirrer was used 25 at low speed (1 to 5 rounds per second) to maintain the cells in an homogenous suspension. The magnetic stirring bar was equipped with a ring to avoid crushing the cells, and the subsequent release of aequorin in the culture medium. Alternatively, a culture vessel equipped for culture of cells in suspension may be used.

One uses the EG&G Wallac's MicroLumat-Plus microplate luminometer, which allows injection and direct subsequent recording of the emitted light from each well of

a 96-well plate. The end of the entrance tube of the dispenser was placed at the bottom of the cell suspension and the dispenser was washed with 3 times the dead volume of the apparatus so that all the volume of the tube and 5 pumps was filled with the cell suspension. The 96-well plate containing the solutions of agonists was then inserted into the luminometer. Then, for each well, 50 μ l of the cell suspension (i.e. 100 000 cells) was dispensed into the well (at the lowest injection speed (0.4 s) to 10 prevent cell breakage that would release aequorin into the culture medium) and the emitted light was immediately recorded during 30 seconds. After reading the first well, cells were injected into the next well and emitted light was recorded, etc. For each plate, a series of curves 15 representing the intensity of the emitted light as a function of time for each well was displayed (figure 1). The intensity of the emitted light was integrated over 30 s using the Winglow software provided with the luminometer, yielding for each well one value representative of the 20 emitted light and hence of the stimulation of the CCR-5 receptor by the agonist present in the well. These values can be plotted against the logarithm of the concentration to generate dose-response curve as shown in figure 2. These allow the determination of half-maximal 25 response doses (EC50) for each ligand. For the generation of these data, 288 measurements were performed in less than 3 hours.

Example 2

A CHO cell line expressing the serotonin 5HT- 2B receptor, the G α 16 coupling protein and apoaequorin was established. Cells were treated as described in Example 1 and were dispensed after dilution (100 μ l/well,

13

corresponding to 50 000 cells) on 100 μ l of solutions of known agonists for this receptor. The emitted light was recorded during 20 s for each well. Dose-response curves obtained for different agonists are shown in figure 3. For the generation of these data, 144 measurements were performed in less than 1 hour.

Example 3

K562 cells expressing the chemokine CCR3

10 receptor were transfected by a plasmid for the expression of aequorin and the Gα16 coupling protein. Cells stably transfected were selected for 2 weeks with the antibiotic Zeocin. These cells were cultured in suspension in DMEM culture medium containing 10% FBS. They were centrifuged and the pellet was used as described in example 1 to perform aequorin measurements. A dose-response curve with eotaxin and MCP-4 generated by this method is described in figure 4.

20 Example 4

orexin 2 receptor, the G α 16 coupling protein and apoaequorin was established. Cells were treated as described in Example 1 and were dispensed after dilution (100 μ 1/well, corresponding to 25 000 cells) on 50 μ 1 of solutions of known agonists for these receptor. The emitted light was recorded during 20 s for each well. A doseresponse curve obtained with this method is showed in figure 5 A.

30

Example 5

A CHO cell line expressing the canabinoid CB1 receptor, the $G\alpha16$ coupling protein and apoaequorin was

14

established. Cells were treated as described in Example 1 and were dispensed after dilution (100 μl/well, corresponding to 50 000 cells) on 100 μl of solutions of a known agonist for this receptor. The emitted light was recorded during 20 s for each well. A dose-response curve obtained with this method is showed in figure 5 B.

For the measure of antagonistic activities, cells were injected on series of dilutions of antagonists. At this stage, the emitted light was recorded to check that 10 the potential antagonists had no agonistic activity. The cells were incubated with the antagonists for 30 min. A solution of agonist (alpha-methyl-5HT) of the receptor was then injected on the cells and the emitted light was immediately recorded for each well. The emitted light was plotted as a function of the logarithm of the 15 antagonist concentration to yield the graph Increasing antagonist concentrations result in a decreasing light emission upon agonist addition. An agonist of another receptor expressed by the cells (usually ATP, acting at P2 receptors) can then be injected on the 20 mixture antagonist, cells and agonists as a control to check that the cells still have active aequorin up to that moment of the experiment. For example, cytotoxic compounds that increase the intracellular calcium concentration can make 25 the aequorin consume all the coelenterazine present in the cell. Such a cytotoxic compound will be detected by the absence of signal upon ATP injection (Fig. 6).

REFERENCES

15

- Button, D. & Brownstein, M., Cell Calcium 14, pp. 663-671 (1993)
- Komatsuzaki, K. et al., FEBS Lett. 406, pp. 165-170
 5 (1997)
 - Milligan, G. et al., TIPS 17, pp. 235-237 (1996)
 - Park, D. et al, J. Biol. Chem. 267, pp. 16048-16055
 (1992).
 - Rink, T.J., FEBS Lett. 268, pp. 381-385 (1990)
- - Sheu, Y.-A. et al, Analytical Biochemistry 209, pp. 343-347 (1993)
- Stables, J. et al., Analytical Biochemistry 252, pp. 115-126 (1997).
 - Tsunoda, Y., Biochim. Biophys. Acta 1154, pp. 105-156 (1993)

CLAIMS

- Diagnostic and/or dosage method of an agonist and/or an antagonist or any modulator for a calcium-coupled receptor or a calcium-coupled channel or any other calcium-coupled protein, comprising the following successive steps:
 - disposing the agonist and/or the antagonist upon a solid support,
- incubating one or more cell(s) expressing apoaequorin or any other related protein and said calcium-coupled receptor with coelenterazine or any other cofactor of a calcium-sensitive protein in order to reconstitute an active aequorin by said cell(s),
- adding to said solid support one or more of said cells,
 and
 - obtaining the measurement of an emitted light by said cell(s).
 - 2. Method according to claim 1, wherein the solid support is a microtiter plate.
- 3. Method according to claim 2, characterised in that the microtiter plate is a 96-well microtiter plate, or a 384-well plate, or a 1536-well-plate or any other format.
- 4. Method according to any one of the preceding claims, characterised in that the cell expresses apoaequorin in the cytoplasm or in the mitochondria or in any other part of the cell.
- 5. Method according to any one of the preceding claims, wherein the cell expressing a calcium-30 coupled receptor is a cell expressing an endogenous or recombinant G-protein-coupled receptor and/or a cell which expresses proteins intended to ensure a coupling of the


17


analysed receptor (endogenous or overexpressed) to the calcium pathway.


- 6. Method according to claim 5, wherein said protein is selected from the group consisting of natural 5 Gα16 or Gα15 protein, chimeric G-protein resulting from a fusion between two different G-proteins or phospholipase Cβ2 protein or any other coupling protein or chemical.
- 7. Method according to any one of the preceding claims, characterised in that the measurement of the emitted light is obtained with one or more luminometer(s), advantageously equipped with several dispensers and measurement heads.
- 8. High-throughput screening diagnostic and/or dosage device intended for the high-throughput screening diagnostic and/or dosage method according to any one of the preceding claims, comprising the following elements:
 - a microtiter plate, preferably a 96-well microtiter plate,
- 20 a medium containing cell(s) expressing apoaequorin and a calcium-coupled receptor,
 - a medium containing coelenterazine, and
 - means for detecting an emitted light by said cell(s).
- 9. Device according to claim 8, comprising 25 means for the automatic performance of the successive steps of the diagnostic and/or dosage method according to any one of the claims 1 to 7.
- 10. Agonist or antagonist of a receptor identified by the method according to any one of the claims 30 1 to 8.

1/6

FIG.1







4/6




FIG. 4

A.

Dose response curve Mix Aequorin-Orexin 1 R

PCT/BE99/00087




FIG. 5A

B.

Dose response curve Mix Aequorin-Orexin 2 R





FIG. 5B



- Methysergide
- Methiothepin
- mesulergine
- Mianserin
- Ketanserine

FIG. 6